

Gait Analysis using a High-Resolution Spatiotemporal Pressure Sensing Insole

Oren Wei, Yao Zhao, Yanhua Chen, Andrea Cheng, Connie Chen, Andrew Ying

Biomedical Pls: Dr. Nitish Thakor, Dr. Ryan Roemmich, Dr. Preeti Raghavan Engineering Pls: Dr. Casey Overby Taylor, Dr. Joseph Greenstein

Department of Biomedical Engineering, Department of Applied Mathematics and Statistics, Whiting School of Engineering

Johns Hopkins University, Baltimore, MD



JOHNS HOPKINS
WHITING SCHOOL
of ENGINEERING


Background

- Gait patterns are **valuable biomarkers** for neurological conditions such as stroke
- Instrumented gait analysis** (motion capture, force plates) is infeasible outside the clinic
- The Thakor Lab has developed an insole sensor with sufficient resolution and sampling rates for **accurate assessment of day-to-day patient gait**

Objectives

- Identify clinically relevant foot regions and derive **novel pressure and gait metrics** from insole data
- Assess metric **interpretability** and **ease-of-adoption** based on literature review and clinician input
- Compare metrics **across walking conditions** for **right leg paretic stroke patient** and **healthy subject**
- Analyze **left-right asymmetries**, **inter-subject variability**, and rank each metric's predictive value

Data Collection Methods



Device: High-resolution (64×16, 100Hz) pressure data collected via sensorized insoles.

Protocol: Split treadmill increasing the speed of the right foot (0.6, 1, 1.4, 1.8 m/s)

Figure 1. Hardware components of the insoles.

Calculating Gait Metrics




Figure 2. Trajectory of Center of Pressure (CoP) over 1 gait cycle.




Figure 3. Peak Pressure (PP) and Regional Impulse (RI) of 2 gait cycles on plot of mean pressure.


- One-Way ANOVA Test conducted for statistical significance between different regions.

Two Methods for Region of Interest Segmentation

Bottom-Up Pressure Valley Segmentation

- Divides foot regions by detecting low-pressure valleys in column-wise and row-wise means
- Smooth (Savitzky-Golay and median filter) and remove outliers

Figure 4a. Pressure valley segmentation. **Figure 4b.** Smoothing filter effects on peak/valley detection.



Clinically Informed Anatomical Segmentation

- Anterior-Posterior axis: Forefoot, Midfoot, Hindfoot (based on joint landmarks)
- Medial and Lateral regions




Figure Citation: Murphy A, Regions of the foot (illustrations). Case study. Radiopaedia.org (Accessed on 20 Apr 2025) <https://doi.org/10.5334/rtd-99469>

Figure 5a. 6 anatomical regions outlined. **Figure 5b.** Locations of joint landmarks.

Comparison of Stroke and Healthy Participants' Gait Using Anatomical Segmentation




Figure 6. PP, RI, and RRI (computed every 10 cycles) for stroke (0.75 m/s) and healthy (0.6 m/s) participants.





Figure 7. CoP trajectories for one gait cycle of healthy and stroke patient. **Figure 8.** CoP instantaneous velocity in anterior-posterior direction average per gait cycle. **Figure 9.** CoP range of right foot in medial-lateral direction.

Key Takeaways

- Figure 6:** stroke patient has significantly greater pressure in all regions except hindlat and hindmed for RRI. **medial-lateral asymmetry** in stroke patient ↑
- Figure 7:** stroke patient **right foot drags behind** left foot
- Figure 8:** stroke patient has reduced CoP velocity at heelstrike and toe off → forward **speed and force upon ground contact is reduced** in stroke
- Figure 9:** Right mediolateral CoP(x) range ↑ → **impaired lateral stability**

Analysis on Healthy Participant Using Pressure Valley Segmentation




Figure 10. Left-right asymmetry in PP/RI/RRI during normal gait (0.6 m/s both feet).




Figure 11. Effects of changing right foot speed on left foot metrics.

Key Takeaways

- Figure 10:** Left foot pressure significantly higher than right across all regions except midfoot for relative region impulse → **left foot acts as the main supporting leg** during walking.
- Figure 11:** Right foot speed ↑ causes **pressure shifts from toes/heels to midfoot in left foot** → enhance stability with balance and propulsion at higher speeds

Conclusions and Next Steps

- High-resolution, high-frequency foot pressure data reveals **meaningful gait differences** across and within individuals
 - These insights **align with clinical expectations** and may support personalized diagnosis and rehabilitation planning
 - This approach shows **promise for precision gait therapy**
- Next Steps:**
- Use dynamic threshold and more conservative filtering
 - Recalculate pressure timing with speed-dependent shifts in mean pressure changes